

Raman spectroscopy as a complementary tool for chitin and chitosan quantification

Dr. Javier Pinto*

Study, Preservation, and Recovery of Archaeological, Historical and Environmental Heritage (AHMat) Research Group, Condensed Matter Physics, Crystallography, and Mineralogy Department, Faculty of Science, University of Valladolid, 47011 Valladolid, Spain
BioEcoUVA Research Institute on Bioeconomy, University of Valladolid, Spain

*Correspondence: javier.pinto@uva.es

Raman spectroscopy has emerged as a powerful technique for biomolecular analysis due to its high chemical and structural sensitivity, non-destructive nature, and lack of sample preparation requirements. This talk will introduce the basic principles of Raman spectroscopy and provide a brief overview of its use in studying fungal cell wall components such as chitin and chitosan. Traditional methods face challenges in distinguishing these compounds, but Raman spectroscopy offers a promising alternative. Recent results will be presented demonstrating the use of Raman spectroscopy to quantify the chitin:chitosan ratio in fungal species such as *Pochonia chlamydosporia* and *Akanthomyces lecanii*.